ATLAS LATERAL MASS SCREW ENTRY POINT USING THE NOVEL INTERSECTION BETWEEN LATERAL MASS AND INFEROMEDIAL EDGE OF THE ATLAS POSTERIOR ARCH TECHNIQUE

Wongthawat Liawrungrueang

School of Medicine, University of Phayao, Phayao, Thailand

Background:

Atlas lateral mass screw placement is a very challenge for Atlas osteosynthesis or combined with atlantoaxial fixation.

Objectives:

To compare the Atlas lateral mass screw placement using the intersection between the lateral mass and inferomedial edge of the posterior arch, as an isolated medial reference, between screw trajectory 0-degree and 15-degree angulation when performing Atlas lateral mass screw fixation.

Materials and methods:

Twenty-four Atlas in fresh cadavers were prepared for inserting the 4.0-mm lateral mass screws. They were divided into two groups: Group 1, screws inserted at 3 mm lateral to the reference entry point with trajectory 0-degree angulation, and Group 2, those inserted with trajectory 15-degree angulation. CT scan was performed before and after inserting Atlas lateral mass screws for evaluating atlas anatomy and screw breach grading.

Results:

All parameters between Group 1 and Group 2 were found statistically different by bilateral intraosseous screw lengths. Twenty-two screws (91.67%) using the 0-degree medial angulation and nineteen screws (79.17%) using the 15-degree medial angulation had no cortical violation (Grade 0). However, two screws (8.33%) with 0-degree medial angulation and five screws (20.83%) with 15-degree medial angulation had breach less than 2 mm (Grade 1) without screw pullout. There was no screw breach distance between 2-4 mm (Grade 2) or breach distance more than 4 mm (Grade 3).

Conclusion:

Bilateral lateral mass screw length in Group 2 was longer than Group 1. Using 3-mm lateral to an intersection between the lateral mass and inferomedial border of the Atlas posterior arch can be used to be an easily-identifiable entry point to insert the lateral mass screw with trajectory 0-degree and 15-degree medial angulation in Atlas fixation.